



Mini Pillow

Fits in the palm of your hand.



FITS IN THE PALM OF YOUR HAND

Mini Pillow is a hand-operated irrigation device that is designed to enhance visibility by efficiently clearing debris during ureterorenoscopy.



DESIGNED FOR THE OPERATOR TO RETAIN CONTROL

- The mini pillow's flat ergonomic design of the irrigation device allows the surgeon to take charge of the irrigation process, as it can be comfortably held against the handle of the ureteroscope for precise control.

LOW VOLUME

- Offers low volume displacement of irrigation fluid preventing large spikes in intrarenal pressure, as seen with higher volume devices.
- A highly responsive device with up to 3ml bolus size, provides instantaneous clearance of the visual field.
- An accelerated, on demand flow achieved with compression of the bulb.

IN FLOW SWITCH

- The In Flow Switch creates a comfortable single-handed on/off control of the flow and prevents inadvertent spillage of irrigation fluid, at the point of connection.

TUBING EXTENSION

- The extension allows the assistant to control the irrigation in a comfortable position.

DEHP FREE

- Makes use of super soft thermoplastic resins, without the plasticizers that harm the environment, and a better choice for the patient too.

Product Code	Description	Qty
URO 052	Mini Pillow	10



Wismed CC
71 Princess Alice Avenue,
Durban, 4001, South Africa
+27 31 206 0177
orders@wismed.co.za
www.wismed.co.za



Wismed PL Sp. z o.o.
ul. Wielka 31/2,
53-338, Wrocław, Poland
+48 882 592 872
biuro@wismedpl.com
www.wismedpl.com



TNMC Devices (Llc)
Ulica Franje Krezme 1
10000 Zagreb, Croatia
Tel: +38 59 1622 9658
Email: tnmcchr@tnmc.biz

