



Smart Bulb

Softer and safer.



SOFTER, SAFER, SMART BULB.

Smart Bulb is a hand held operated irrigation device that provides cooling and enhanced visibility by efficiently cleaning debris during ureterorenoscopy.

ADVANTAGES

SOFTNESS + WIDE OUTFLOW TUBE = REDUCED FATIGUE

- Softness of the material, together with a wider outflow tube, reduces the flow resistance during hand held controlled irrigation, which minimizes operator's fatigue.
- An accelerated, on demand flow achieved with compression of the bulb effortlessly.

IN FLOW SWITCH

- This user-friendly In Flow Switch creates a comfortable single-handed on/off control of the flow and prevents inadvertent spillage of irrigation fluid at the point of connection.

TUBING EXTENSION

- The extension allows the assistant to control the irrigation in a comfortable position.

DEHP FREE

- Makes use of super soft thermoplastic resins, without the plasticizers that harm the environment, and a better choice for the patient too.

| Product Code | Description | Qty |
|--------------|-------------|-----|
| URO 037 | Smart Bulb | 10 |



Wismed CC
71 Princess Alice Avenue,
Durban, 4001, South Africa
+27 31 206 0177
orders@wismed.co.za
www.wismed.co.za



Wismed PL Sp. z o.o.
ul. Wielka 31/2,
53-338, Wrocław, Poland
+48 882 592 872
biuro@wismedpl.com
www.wismedpl.com



TNMC Devices (Llc)
Ulica Franje Krezme 1
10000 Zagreb, Croatia
Tel: +38 59 1622 9658
Email: tnmcchr@tnmc.biz

